

EXHIBIT B

LABORATORY NOTEBOOK

GENENTECH, INC.

26865



TITLE _____

From Page No. _____

REDACTEDSet up TX of 293S cells. 100mm Dishes.
Feed cells 8:00am.Tube A
0.5 ml's 1:10 TE
50 µl 2.5M CaCl₂
1 µl VAPNA
10 µg DNA BelowTube B
0.5 ml's 2x HRSDNA on cells 10:00am
Feed Fresh media 9:00pm

Plate # DNAs

1	PrKS	
2	HVEM	
3	TNFR1	
4	Apo-3	
5	PrKS	2 µg + Ap-1
6	HVEM	Luciferase
7	TNFR1	
8	APO-3	
9	PrKS	2 µg + NFKB
10	HVEM	Luciferase
11	TNFR1	
12	APO-3	
13	PrKS	
14	DO 1.1	Enzyme
15	" 2.1	2098364
16	" 3.1	1237537
17	" 4.1	

For C-Jun (SAPK Assay)

Go To pg # 51

For Luciferase
Below.Gave To Tamie
For Facs.

2.1 kills cells.

Harvest Luciferase cells 5-12 in 1ml
cell lysis Buffer. put on rotator 30 min.
Count on Luciferase Assay.
To pg # 49

To Page No. _____

Witnessed & Understood by me, _____

Date _____

Invented by _____

Date _____

Recorded by _____